



Leica MacroFluoTM

Technical Information

Living up to Life

Leica
MICROSYSTEMS





Leica MacroFluo™ System, equipped with Leica Z16 APO A macroscope, Leica DFC camera, motorized zoom, transmitted light base TL RCI™, SmartTouch™ control panel, and Leica EL6000 fluorescent light source.

Leica Z6 APO & Z16 APO / Z6 APO A & Z16 APO A

Technical data, performance characteristics

Designation	Leica MacroFluo™ with 6.3:1 zoom	Leica MacroFluo™ with 16 : 1 zoom
Type of microscope	Macroscope with vertical beam path, apochromatic zoom system, incident illuminator, fluorescence illumination coupled with the optical zoom system, and trinocular tube	
Zoom system	Apochromatic 6.3 : 1 zoom Leica Z6 APO, lead-free	Apochromatic 16 : 1 zoom Leica Z16 APO, lead-free
Zoom factor	0.57 – 3.6×	0.57 – 9.2×
Engageable zoom positions	0.57 / 0.8 / 1 / 1.25 / 1.6 / 2 / 2.5 / 3.2 / 3.6	0.57 / 0.8 / 1 / 1.25 / 1.6 / 2 / 2.5 / 3.2 / 4 / 5 / 6.3 / 8 / 9.2
Tubes	Binocular phototube HC L2TU 4/5/7 with photo-TV output 100%:0% / 0%:100%, 20° viewing angle Trinocular tube BDT 25+V100/50/0 with photo-TV output 100%:0% / 50%:50% / 0%:100%, 30° viewing angle Tube BDTP 25 100/50/0 with fixed port, photo-TV output 100%:0% / 50%:50% / 0%:100%, 30° viewing angle	
Built-in iris diaphragm	For continuous depth of field control	
Designation	Leica MacroFluo™ with automated 6.3 : 1 zoom	Leica MacroFluo™ with automated 16 : 1 zoom
Type of microscope	Macroscope with vertical beam path, motorized, apochromatic zoom system, incident illuminator, fluorescence illumination coupled with the optical zoom system, and trinocular tube	
Zoom system, motorized	Apochromatic 6.3 : 1 zoom Leica Z6 APO A, lead-free	Apochromatic 16 : 1 zoom Leica Z16 APO A, lead-free
Zoom factor	0.57 – 3.6×	0.57 – 9.2×
Motorized zoom	Continuously motorized and approachable positions	Continuously motorized and approachable positions
Tubes	Binocular phototube HC L2TU 4/5/7 with photo-TV output 100%:0% / 0%:100%, 20° viewing angle Trinocular tube BDT 25+V100/50/0 with photo-TV output 100%:0% / 50%:50% / 0%:100%, 30° viewing angle Tube BDTP 25 100/50/0 with fixed port, photo-TV output 100%:0% / 50%:50% / 0%:100%, 30° viewing angle	
Motorized iris diaphragm	For continuous depth of field control	For continuous depth of field control
Motorized fine focus	Integrated, 10 mm travel	Integrated, 10 mm travel
Fluorescence illuminator		
Type	Coaxial incident illuminator L4/23 and carrier for incident light fluorescence and incident light brightfield with: • Reflector disc for 5 filter blocks, 0-pixel technology = no pixel offset when changing • Aperture and field diaphragm can be centered • Switchable blue filter BG 38 and UV light stop • Light trap for suppression of ambient light	
Fluorescence filter	• Filter system (excitation filter, beam-splitting mirror, and barrier filter) on a filter block • Extensive filter block selection, filter combinations based on customer specifications	
Light source	Leica EL6000 external fluorescent light source EL6000 - External light source with liquid light guide - Adjustment-free metal halide lamp - Adjustable intensity of the excitation light 106 z lamp housing with: 100 W high-pressure mercury burner (direct current, stabilized/nonstabilized, type 103 W/2) - Chromatically-corrected collector, focusable - Lamp mount can be centered - Supply unit ebq 100	
Safety measures	UV protection screen, UV light stop, stray light protection for lamp housing	

Optical accessories			
Objectives	Type/magnification factor	Working distance	Max. object field Ø
	• Achromatic 0.32×	297 mm	112 mm
	• Achromatic 0.63×	149 mm	56 mm
	• Planapochromatic 0.5×	187 mm	70 mm
	• Planapochromatic 0.8×	112 mm	44 mm
	• Planapochromatic 1×	97 mm	35 mm
	• Planapochromatic 2×	39 mm	17.5 mm
	• Planapochromatic 5×/ 0.5 LWD	19 mm	5 mm
Objective adapters	• For Leica M-series achromatic objectives • For 10× and 20× Leica DM objectives		
DM objectives (only for upper zoom range)	• DM objective N Plan L 20× / 0.40 corr • DM objective N Plan 10× / 0.25– / A5.8		
Fine focusing	10 mm travel, 1 micron resolution (Leica Z6 APO / Leica Z16 APO) (motorized and integrated for Leica Z6 APO A and Leica Z16 APO A)		
Interpupillary distance	55 mm – 75 mm		
Wide-field eyepieces for spectacle wearers	Widefield HC plan eyepieces 10×/20, 10×/22, 10×/25		
Imaging			
Image acquisition, storage, and processing	Various digital FireWire camera systems for fluorescence		
Software for image archiving, analysis, and processing	LAS (Leica Application Suite), LAS AF (Leica Z6 APO A and Leica Z16 APO A motorized zooms), QWin, Leica MM AF powered by Metamorph®		
Tubes	• Binocular phototube HC L2TU 4/5/7 • Tube BD TP 25 100/50/0 with fixed port • Tube BD T 25+ V100/50/0, mechanical		
Stands, illuminators			
Transmitted light stands	Transmitted light bases TL ST, TL BFDF, TL RC™, and TL RCI™ for a variety of applications and illumination options		
Incident light stand	With stage plate, black and white, 120 mm		
Motorized focus system	With 620 mm column, control via handswitch, SmartTouch™, footswitch, or computer		
Manual focusing drive	Coarse/fine focus, adjustable ease of movement, with 620 mm column		
Stages	Leica IsoPro™ Motorized and Manual Cross-stage, gliding stage, Leica MATS Thermocontrol System heating stage, polarization		
Additional illuminators	Various LED illuminators and cold light sources		

For additional accessories, see assembly diagram and scope of delivery.

Optical data – with HC plan eyepieces 10×/22

Objectives	Planapochromatic 1×		Planapochromatic 0.5×		Planapochromatic 0.8×		Planapochromatic 2×		Planapochromatic 5×		
Working distances	97 mm		187 mm		112 mm		39 mm		20 mm		
Eyepieces	Total magnification	Object field ∅	Total magnification	Object field ∅	Total magnification	Object field ∅	Total magnification	Object field ∅	Total magnification	Object field ∅	
Zoom position	×	mm	×	mm	×	mm	×	mm	×	mm	
6.3:1	0.57	7.1	30.9	3.6	61.8	5.7	38.6	14.3	15.4	35.6	6.18
	0.8	10.0	22.0	5.0	44.0	8.0	27.5	20.0	11.0	50.0	4.40
	1	12.5	17.6	6.3	35.2	10.0	22.0	25.0	8.80	62.5	3.52
	1.6	20.0	11.0	10.0	22.0	16.0	13.8	40.0	5.50	100	2.20
	2	25.0	8.80	12.5	17.6	20.0	11.0	50.0	4.40	125	1.76
	2.5	31.3	7.04	15.6	14.1	25.0	8.80	62.5	3.52	156	1.41
	3.2	40.0	5.50	20.0	11.0	32.0	6.88	80.0	2.75	200	1.10
	3.6	45.0	4.89	22.5	9.78	36.0	6.11	90.0	2.44	225	0.98
	4	50.0	4.40	25.0	8.80	40.0	5.50	100	2.20	250	0.88
	5	62.5	3.52	31.3	7.04	50.0	4.40	125	1.76	313	0.70
	6.3	78.8	2.79	39.4	5.59	63.0	3.49	158	1.40	394	0.56
8	100	2.20	50.0	4.40	80.0	2.75	200	1.10	500	0.44	
9.2 (16:1)	115	1.91	57.5	3.83	92.0	2.39	230	0.96	575	0.38	

HC plan eyepieces 10×/22, binocular phototube HC L2TU 4/5/7 (factor 1.25×), open iris diaphragm

Data with digital camera Leica DFC310 FX Full Frame

Objectives	Planapochromatic 1×			Planapochromatic 0.5×			Planapochromatic 0.8×			Planapochromatic 2×			Planapochromatic 5×		
Working distances	97 mm			187 mm			112 mm			39 mm			20 mm		
Eyepieces	Magnification chip: specimen	Object field projected onto chip		Magnification chip: specimen	Object field projected onto chip		Magnification chip: specimen	Object field projected onto chip		Magnification chip: specimen	Object field projected onto chip		Magnification chip: specimen	Object field projected onto chip	
Zoom position		mm × mm			mm × mm			mm × mm			mm × mm			mm × mm	
6.3:1	0.57	0.50	18.0 13.4	0.25	36.0 26.9	0.40	22.5 16.8	1.00	9.00 6.72	2.49	3.60 2.69				
	0.8	0.70	12.8 9.58	0.35	25.7 19.2	0.56	16.0 12.0	1.40	6.41 4.79	3.50	2.57 1.92				
	1	0.88	10.3 7.67	0.44	20.5 15.3	0.70	12.8 9.58	1.75	5.13 3.83	4.38	2.05 1.53				
	1.6	1.40	6.41 4.79	0.70	12.8 9.58	1.12	8.02 5.99	2.80	3.21 2.40	7.00	1.28 0.96				
	2	1.75	5.13 3.83	0.88	10.3 7.67	1.40	6.41 4.79	3.50	2.57 1.92	8.75	1.03 0.77				
	2.5	2.19	4.10 3.07	1.09	8.21 6.13	1.75	5.13 3.83	4.38	2.05 1.53	10.9	0.82 0.61				
	3.2	2.80	3.21 2.40	1.40	6.41 4.79	2.24	4.01 2.99	5.60	1.60 1.20	14.0	0.64 0.48				
	3.6	3.15	2.85 2.13	1.58	5.70 4.26	2.52	3.56 2.66	6.30	1.43 1.06	15.8	0.57 0.43				
	4	3.50	2.57 1.92	1.75	5.13 3.83	2.80	3.21 2.40	7.00	1.28 0.96	17.5	0.51 0.38				
	5	4.38	2.05 1.53	2.19	4.10 3.07	3.50	2.57 1.92	8.75	1.03 0.77	21.9	0.41 0.31				
	6.3	5.51	1.63 1.22	2.76	3.26 2.43	4.41	2.04 1.52	11.0	0.81 0.61	27.6	0.33 0.24				
8	7.00	1.28 0.96	3.50	2.57 1.92	5.60	1.60 1.20	14.0	0.64 0.48	35.0	0.26 0.19					
9.2 (16:1)	8.05	1.12 0.83	4.03	2.23 1.67	6.44	1.39 1.04	16.1	0.56 0.42	40.3	0.22 0.17					

Binocular phototube HC L2TU 4/5/7 (factor 1.25×), video objective 0.7×, open iris diaphragm



Z6 APO / Z6 APO A: Zoom 0.57× – 3.6×
Z16 APO / Z16 APO A: Zoom 0.57× – 9.2×

Leica MacroFluo™ with 6.3:1 zoom (Leica Z6 APO or Leica Z6 APO A)

Data visual with binocular phototube HC L2TU 4/5/7 (tube factor 1.25×) / eyepieces 10×/22			
Objective	Planapochromatic 1×	Planapochromatic 2×	Planapochromatic HR 5×
Zoom range	7.1× – 45×	14.3× – 90×	35.5× – 225×
Resolution	60 – 351 Lp/mm	120 – 702 Lp/mm	max. 1500 Lp/mm
Max. numerical aperture	0.117	0.234	0.5
Finest visible structure width	1.4 µm	0.71 µm	0.33 µm
Object field ∅	4.9 – 30.9 mm	2.4 – 15.4 mm	1.0 – 6.2 mm
Max. depth of field	3.1 mm	0.8 mm	0.1 mm
Working distance	97 mm	39 mm	19 mm
Data with digital camera Leica DFC310 FX Full Frame / binocular phototube HC L2TU 4/5/7 (tube factor 1.25×) / video objective 0.7×			
Magnification chip: Object	0.5× – 3.15×	1.0× – 6.3×	2.5× – 15.8×
Digital resolution	19.3 – 122 Lp/mm	35 – 220 Lp/mm	97 – 610 Lp/mm
Object field projected onto chip	18 × 13.4 mm / 2.9 × 2.1 mm	9 × 6.7 mm / 1.43 × 1.06 mm	3.6 × 2.7 mm / 0.57 × 0.4 mm
Max. depth of field	1.4 mm	0.33 mm	0.05 mm

Leica MacroFluo™ with 16:1 zoom (Leica Z16 APO or Leica Z16 APO A)

Data visual with binocular phototube HC L2TU 4/5/7 (tube factor 1.25×), eyepieces 10×/22			
Objective	Planapochromatic 1×	Planapochromatic 2×	Planapochromatic HR 5×
Zoom range	7.1× – 115×	14.3× – 230×	35.5 – 575×
Resolution	51 – 336 Lp/mm	102 – 672 Lp/mm	max. 1500 Lp/mm
Max. numerical aperture	0.112	0.224	0.5
Finest visible structure width	1.49 µm	0.75 µm	0.33 µm
Object field ∅	1.9 – 30.9 mm	0.96 – 15.4 mm	0.38 – 6.2 mm
Max. depth of field	3.8 mm	0.94 mm	0.15 mm
Working distance	97 mm	39 mm	19 mm
Data with digital camera Leica DFC310 FX Full Frame / binocular phototube HC L2TU 4/5/7 (tube factor 1.25×) / video objective 0.7×			
Magnification chip: Object	0.5× – 8.05×	1.0× – 16.1×	2.5× – 40.3×
Digital resolution	19.3 – 312 Lp/mm	38.7 – 624 Lp/mm	96.7 – 1500 Lp/mm
Object field projected onto chip	18 × 13.4 mm / 1.1 × 0.8 mm	9 × 6.72 mm / 0.56 × 0.42 mm	3.6 × 2.7 mm / 0.2 × 0.17 mm
Max. depth of field	1.7 mm	0.35 mm	0.07 mm



Name	MPixel	Sensor	Size	Pixel-Size	Cooling	Interface	C-Mount
DFC295	3	CMOS-c	1/2"	3.2 µm	no	FWb	0.5x
DFC290 HD	3	CMOS-c	1/2"	3.2 µm	no	FWb+HDMI	0.5x
DFC400	1.4	CCD-c	1/2"	4.65 µm	no	FWb+Sync	0.5x
DFC450	5	CCD-c	2/3"	3.4 µm	no	FWb	0.63x
DFC450 C	5	CCD-c	2/3"	3.4 µm	yes	FWb+Sync	0.63x
DFC495	8	CCD-c	2/3"	2.7 µm	yes	FWb+Sync	0.63x
DFC310 FX	1.4	CCD-c	2/3"	6.45 µm	yes	FWb+Sync	0.63x
DFC345 FX	2	CCD-m	1/2"	4.4 µm	yes	FWb+Sync	0.5x
DFC365 FX	1.4	CCD-m	2/3"	6.45 µm	yes	FWb+Sync	0.63x
DFC500	12	CCD-c	2/3"	6.45 µm	yes	FWa+Sync	0.63x

-c: color / -m: monochrome

V2 – 2012

Digital camera systems

- 12 730 469 Leica DFC295 camera kit
- 12 730 468 Leica DFC290 HD camera kit*
- 12 730 470 Leica DFC400 camera kit*
- 12 730 411 Leica DFC450 camera kit*
- 12 730 412 Leica DFC450 C camera kit*
- 12 730 471 Leica DFC495 camera kit*
- 11 547 001 Leica DFC310 FX camera kit for Mac
- 11 547 002 Leica DFC310 FX camera kit for PC
- 11 547 003 Leica DFC345 FX camera kit*
- 11 547 004 Leica DFC365 FX camera kit
- 12 730 066 Leica DFC500 camera kit

* Mac-Support

Digital camera system accessories

- 12 446 120 FireWire cable, FW-a-a, 2.5 m, 6-pin to 6-pin
- 12 730 186 FireWire cable, FW-b-b, 2.5 m, 9-pin to 9-pin
- 12 730 187 FireWire cable, FW-a-b, 2.5 m, 6-pin to 9-pin
- 12 730 211 HDMI cable, 3m, HDMI both ends
- 12 730 183 OHCI FireWire PCI32 card, FWa+b+USB2
- 12 730 210 OHCI FireWire PCI express card, FW-b
- 12 730 446 FW-b Notebook Kit (includes PCcard Express, power supply: 100–240 V, 24 W, adapter FW-b-a)
- 12 730 447 OHCI FireWire PCI express card Low Profile, FW-b
- 12 730 225 PowerKit for DFC290 HD, FW-b

Phototubes and camera adapter


- 11 621 500 MPV adapter, HC/delta
- 11 541 038 Vario adapter for C-mount 0.55x to 1.1x – Delta
- 11 541 517 Vario adapter for C-mount 0.32x to 1.6x – HC
- 11 543 702 B-mount adapter 1x
- 11 543 706 C-mount adapter 1x
- 11 541 543 C-mount adapter HC 0.7x, for 2/3" digital cameras
- 11 541 544 C-mount adapter HC 0.55x, for 1/2" digital cameras
- 11 541 510 C-mount adapter HC 1x, for 1" digital cameras

For connection to the camera diagram, see previous page


-
- 10x**
11 507 802
11 507 808
11 507 807
- 10x**
11 501 598
11 505 162
11 505 146
(for 11 505 146)
- Filtercube k-size**
11 505 235
- 11 504 117**
- 11 504 116**
- 11 504 115**
- RS232 – USB**
11 532 295
- 11 504 101**
- 10 399 211**
10 446 421
- 10 446 154**
- Z6 APO**
10 447 174
- Z16 APO**
10 447 173
- Z6 APO A**
10 446 368
- Z16 APO A**
10 446 369
- M12 cable, 10 m**
10 447 469
- 4**
- Distribution box**
For Leica Z-series
- Foot switch**
10 447 398
- RS232**
10 447 401
- CTL2 cable, 2 m**
10 447 379
- Smartfouch™**
10 450 052
- 2**
- Objective slide**
For Leica Z6 APO &
amp; Leica Z16 APO
10 450 094
10 450 138
- Fine focus**
10 447 175

Objectives and optical accessories

10 447 177	Planapochromatic objective 0.5×
10 446 360	Planapochromatic objective 0.8×
10 447 176	Planapochromatic objective 1×
10 447 178	Planapochromatic objective 2×
10 447 243	Planapochromatic objective 5×
10 450 332	Adapter for achromatic objectives for M series
10 450 191	Achromatic objective 0.32×
10 450 192	Achromatic objective 0.5×
10 450 160	Achromatic objective 0.63×
10 450 243	Adapter for LWD achromatic objectives
10 411 597	Achromatic objective f = 100 mm
10 441 787	Achromatic objective f = 150 mm
10 431 692	Achromatic objective f = 175 mm
10 382 162	Achromatic objective f = 200 mm
10 382 172	Achromatic objective f = 400 mm
35 000 011	Achromatic objective f = 800 mm
35 000 013	Achromatic objective f = 1200 mm
10 447 179	Adapter Leica DM objectives to Leica Z6/Leica Z16
10 450 169	Leica LED3000 NVI™ – vertical illuminator for using objectives of 58 mm maximum diameter
10 450 271	LED3000 RL – ring illuminator, Ø 58 mm objectives, 24 power LEDs, 5600 K color temperature, optimized LED auxiliary lens, selectable segments, optimized for 60–150 mm working distance
10 450 548	LED5000 SLI™, spotlight illumination, double-armed gooseneck 500 mm long, with 2 power LEDs, 5600 K color temperature, control unit on separate gooseneck, incl. diffuser pair
10 450 561	Leica LED5000 MCI™ (Multi Contrast Illumination): three adjustable illuminator arcs, each with three LEDs for complete integration with the stereomicroscope; modes fully reproducible with Leica Application Suite including automatic changing of multiple illumination scenarios.
10 450 266	Power supply unit for Leica LED3000/Leica LED5000
35 000 031	Protective glass for planapochromatic objectives
10 315 306	Analyzer in revolving mount
10 447 207	Quarter-wave plate, Ø 58 mm
10 447 460	Quarter-wave plate for 2.0× and 5.0×
10 450 049	Incident-light and transmitted-light bases
10 450 049	Incident-light base
10 450 123	Transmitted-light base TL3000 ST
10 450 124	Transmitted-light base TL4000 BFDF
10 450 125	Transmitted-light base TL4000 RC™ for external cold light sources
10 450 126	Transmitted-light base TL4000 RC™ with integr. halogen illumination
10 450 541	Transmitted-light base TL5000 Ergo with integr. LED illumination
10 450 260	Universal base XL for specimens up to 300 × 300 mm
10 450 297	Adapter for base 10 450 260 for all swinging-arm columns
10 450 562	Stages
10 450 562	Standard stage for transmitted-light bases TL4000 BFDF, TL4000 RC™ and TL4000 RC™




- 10 450 127 Leica IsoPro™ Manual Cross-stage for transmitted-light bases
TL4000 BFDf, TL4000 RC™, TL4000 RCi™ and incident-light base
(with 10 450 122 adapter)
- 10 450 218 Leica IsoPro™ Motorized Cross-stage for transmitted-light bases
TL4000 BFDf, TL4000 RC™, TL4000 RCi™ and incident-light base
(with 10 450 122 adapter)
- 10 450 122 Adapter between cross-stage and incident-light base 10 450 049
- 10 447 275 Leica MATS TL heating stage insert with control unit for TL bases
- 10 447 276 Adapter for stages with Ø 120 mm
- 10 447 391 Stage for LiveOnStage accessories
- 10 447 392 Universal carrier for Petri dishes, specimen slides etc.
- 10 446 301 Gliding stage, Ø 120 mm
- 10 446 302 Polarization stage, Ø 120 mm
- 10 382 130 Attachable mechanical stage for polarization stage
- 10 361 719 Sensitive-tint plate for Pol rotating stage
- 10 446 303 Cup stage, Ø 120 mm
- 33 000 600 Cup stage, Ø 120 mm, surface Ø 150 mm, revolving
- 10 446 228 Glass insert with Pol, Ø 120 mm
- 10 447 057 Stage plate b/w, Ø 120 mm
- 10 447 058 Clear glass insert, Ø 120 mm
- 10 447 153 Matte glass insert, Ø 120 mm
- 10 450 058 Stage plate, b/w for TL bases
- 10 450 059 Additional buttons for IsoPro™ Manual Cross-stage
- Focusing drive**
- 10 450 504 Coarse/fine focusing drive with profile column 420 mm
- 10 450 505 Coarse/fine focusing drive with profile column 620 mm
- 10 450 502 Motorized focus with profile column 420 mm
- 10 450 503 Motorized focus with profile column 620 mm
- 10 450 289 XL extension – for viewing large specimens
- 10 450 205 Combination light guide adapter for installation on the focus column
- 10 447 379 CTL2 cable, 2 m
- Filters**
- 10 447 400 Daylight filter for TL ST base
- 10 447 394 BG38 filter for TL4000 RC™/RCi™ transmitted-light base
- 10 447 395 UV filter for TL4000 RC™/RCi™ base
- 10 447 393 ND filter (gray filter) for TL RC™/RCi™ base
- 10 450 079 Daylight filter for TL4000 RCi™ base
- 10 450 563 Milk glass filter for TL5000 Ergo
- 10 450 564 Green filter for TL5000 Ergo
- 10 450 565 Polarization filter for TL5000 Ergo
- Controls**
- 11 505 180 Leica Smart Move control unit for Leica IsoPro™ motorized cross-stage
- 10 447 444 Leica PSC control unit for Leica IsoPro™ motorized cross-stage
- 10 450 298 Manual controller for motorized focus
- Ergonomic accessories
- 10 447 431 Leica ErgoRest™ handrest for fatigue-free working




Leica Z6 APO & Z16 APO

Dimensions of zoom system, coaxial incident light housing




Leica Z6 APO









Leica Z16 APO

Coaxial incident light housing








Leica Z6 APO A & Z16 APO A


Dimensions of zoom systems




Dimensions of objectives, adapters



Adapter for achromatic objectives



Adapter between fluorescent bridge and focus column



The statement by Ernst Leitz in 1907, **"With the User, For the User,"** describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: **Living up to Life.**

Leica Microsystems operates globally in four divisions, where we rank with the market leaders.

LIFE SCIENCE DIVISION

The Leica Microsystems Life Science Division supports the imaging needs of the scientific community with advanced innovation and technical expertise for the visualization, measurement, and analysis of microstructures. Our strong focus on understanding scientific applications puts Leica Microsystems' customers at the leading edge of science.

INDUSTRY DIVISION

The Leica Microsystems Industry Division's focus is to support customers' pursuit of the highest quality end result. Leica Microsystems provide the best and most innovative imaging systems to see, measure, and analyze the microstructures in routine and research industrial applications, materials science, quality control, forensic science investigation, and educational applications.

BIOSYSTEMS DIVISION

The Leica Microsystems Biosystems Division brings histopathology labs and researchers the highest-quality, most comprehensive product range. From patient to pathologist, the range includes the ideal product for each histology step and high-productivity workflow solutions for the entire lab. With complete histology systems featuring innovative automation and Novocastra™ reagents, Leica Microsystems creates better patient care through rapid turnaround, diagnostic confidence, and close customer collaboration.

MEDICAL DIVISION

The Leica Microsystems Medical Division's focus is to partner with and support surgeons and their care of patients with the highest-quality, most innovative surgical microscope technology today and into the future.

Leica Microsystems – an international company with a strong network of worldwide customer services:

Active worldwide		Tel.	Fax
Australia · North Ryde	+61	2 8870 3500	2 9878 1055
Austria · Vienna	+43	1 486 80 50 0	1 486 80 50 30
Belgium · Groot Bijgaarden	+32	2 790 98 50	2 790 98 68
Canada · Concord/Ontario	+1	800 248 0123	847 405 0164
Denmark · Ballerup	+45	4454 0101	4454 0111
France · Nanterre Cedex	+33	811 000 664	1 56 05 23 23
Germany · Wetzlar	+49	64 41 29 40 00	64 41 29 41 55
Italy · Milan	+39	02 574 861	02 574 03392
Japan · Tokyo	+81	3 5421 2800	3 5421 2896
Korea · Seoul	+82	2 514 65 43	2 514 65 48
Netherlands · Rijswijk	+31	70 4132 100	70 4132 109
People's Rep. of China · Hong Kong	+852	2564 6699	2564 4163
· Shanghai	+86	21 6387 6606	21 6387 6698
Portugal · Lisbon	+351	21 388 9112	21 385 4668
Singapore	+65	6779 7823	6773 0628
Spain · Barcelona	+34	93 494 95 30	93 494 95 32
Sweden · Kista	+46	8 625 45 45	8 625 45 10
Switzerland · Heerbrugg	+41	71 726 34 34	71 726 34 44
United Kingdom · Milton Keynes	+44	800 298 2344	1908 246312
USA · Buffalo Grove/Illinois	+1	800 248 0123	847 405 0164